

Six Transmitter Data communication Manual



If you are using our evaluation software "bioMON" this document is not relevant for you.

You may download "bioMON" as a Windows installer package at www.jobst-technologies.com/download/software/bioMON4_Setup.zip and get advice for installation and operation watching the video tutorials available at www.jobst-technologies.com/support-2/biomon-video-tutorials/.

Six transmitter data communication.....	2
USB version.....	2
OEM version.....	2
Data conversion.....	3
Serial Protocol.....	4
Data Telegram format.....	4
Error Telegram format.....	5
Application example.....	6
Example Pseudo Code on python basis.....	7

Disclaimer: Not for medical, diagnostics and usage on humans. For evaluation use only.

Six transmitter data communication

You will have either the regular version with an USB cable or the OEM version without housing and cable.

USB version

Win 10

Follow these steps:

1. Assure you have Internet access.
2. Plug in the USB cable to the PC.
3. Download our evaluation software "bioMON":
https://www.jobst-technologies.com/download/software/bioMON4_Setup.zip
4. Install bioMON (this requires Admin privileges). Windows 10 will then automatically install a serial-to-USB driver for SIX.
5. You can keep bioMON installed, no matter if you will use it or not.

Win 7

Follow these steps:

1. Download the serial-to-USB driver from:
http://www.prolific.com.tw/US/ShowProduct.aspx?p_id=225&pcid=41
that provides you a virtual COM port.
2. Plug in the USB cable to the PC.

OEM version

For the OEM version connect ground and supply voltage (3.3 – 5.5 V) to your power supply. Connect the Read data (Rx) pin of your microprocessor to the Send data (Tx) pin of the Six transmitter. The pin assignment can be found in the instrumentation datasheet

Biosensor Transmitter OEM:

https://www.jobst-technologies.com/download/datasheets/instrumentation-datasheets/PDD0113-1_Transmitter_OEM_datasheet_en.pdf

There is no need to connect the Rx pin of the Six transmitter as the data is sent by the transmitter automatically every 1.7 seconds.

For the virtual serial port via USB interface or the TTL serial port (UART) with OEM version use port settings of 9600 baud, no parity, 8 data bits, and 1 stop bit (9600,n,8,1).

Disclaimer: Not for medical, diagnostics and usage on humans. For evaluation use only.

Serial Protocol

Data Telegram format

byte#	hex	Description	Comments	
1.	68	Start		
2.	nn	Length	The length byte is the length from the number of data bytes + 1 (the message type id byte). Decimal value is 19.	
3.	nn	Length		
4.	68	Start		
5.	04	Message type ID	Message type id is a unique number identifying the which message should be used to interpret the frame. Value is 4.	
6.	xx	Channel1 Hi Byte	<p>Frame data.</p> <p>If the value of one of the six channels is above 2^{15} counts this represents a negative number. Then subtract 2^{16}.</p> <p>That gives a reading range of (-2^{15}) to $+(2^{15}-1)$ counts. That range is equivalent to either $\pm 25nA$ or $\pm 50nA$ (check your Six transmitters label).</p> <p>If the values are exactly 32767 or -32768 this indicates a value above/below the measurement boundaries i.e. an error.</p>	
7.	xx	Channel1 Lo Byte		
8.	xx	Channel2 Hi Byte		
9.	xx	Channel2 Lo Byte		
10.	xx	Channel3 Hi Byte		
11.	xx	Channel3 Lo Byte		
12.	xx	Channel4 Hi Byte		
13.	xx	Channel4 Lo Byte		
14.	xx	Channel5 Hi Byte		
15.	xx	Channel5 Lo Byte		
16.	xx	Channel6 Hi Byte		
17.	xx	Channel6 Lo Byte		
18.	xx	Temperature Hi Byte		Divide the value by 16 to get the temperature in deg C.
19.	xx	Temperature Lo Byte		
20.	xx	ID Highest Byte (MSB)		
21.	xx	ID Byte 2		
22.	xx	ID Byte 3		
23.	xx	ID Lowest Byte (LSB)		
24.	xx	Checksum	The checksum is the last byte of sum of the values of the bytes of the message type id and the frame data (Bytes 5-23).	
25.	16	Stop byte		


Disclaimer: Not for medical, diagnostics and usage on humans. For evaluation use only.

Error Telegram format

1.	68	Start	
2.	nn	Length	The length byte is the length from the number of data bytes + 1 (the message type id byte). Decimal value is 2.
3.	nn	Length	
4.	68	Start	
5.	05	Message type ID	Value for Error message is 5.
6.	xx	Error code	
7.	xx	Checksum	The checksum is the last byte of sum of the values of the bytes of the message type id and the frame data (Bytes 5-6).
8.	16	Stop byte	

Disclaimer: Not for medical, diagnostics and usage on humans. For evaluation use only.

Application example



$$\text{For mM reading: Glucose_Signal_1 [mM]} = (\text{Chn2}-\text{Chn1}) * \text{G_Glucose1} / 100 / \exp(0.038 * (\text{Temperature} - 32))$$

$$\text{For nA reading: Glucose_Signal_1 [nA]} = (\text{Chn2}-\text{Chn1}) * \text{G_Blank1} / 100 / \exp(0.038 * (\text{Temperature} - 32))$$

#G_Blank1 is the blank gain factor provided in the .def file: $\text{G_Blank1} = 50 / ((2^{15}) - 1) * 100 = 0.1526$

Disclaimer: Not for medical, diagnostics and usage on humans. For evaluation use only.

Example Pseudo Code on python basis

The basic code example receives the byte array and gives the current in nA (for a 50 nA transmitter) per channel as well as the temperature and a timestamp

```
1. # -*- coding: utf-8 -*-
2. For Python3.7
3. @author: Jobst Technologies GmbH
4.
5. This simple script reads data from a COM port, to which a SIX Potentiostat is
6. connected. Then it stores the timestamp in s, the current in nA and the temperature in de-
7. gree C in a tab-delimited txt file.
8. The SIX potentiostat delivers a data point every 1.7 seconds, which is shown on
9. screen and stored a tab-delimited txt file:
10.     Timestamp    in s
11.     ch1 current in nA
12.     ch2 current in nA
13.     ch3 current in nA
14.     ch4 current in nA
15.     ch5 current in nA
16.     ch6 current in nA
17.     Temperature in °C
18.
19. Start conditions:
20. - SIX potentiostat is connected to a known COM port
21. - output file name
22. ""
23. import serial
24. import time
25.
26. # Input data:
27. COM_PORT = "COM5"           # This needs to be changed accordingly
28. OUT_FILENAME = "out_data.txt" # This needs to be changed accordingly
29.
30. # Serial port settings
31. BAUD = 9600
32. TIMEOUT = 0.5
33.
34. # Data telegram settings
35. package_length = 25 # the SIX potentiostat delivers 25-bytes messages
36. data_block = [b'\x00'] * package_length # ring array for incoming data
37. start_timestamp = None
38.
39. # Starting serial port connection
40. with serial.Serial(COM_PORT, baudrate=BAUD, timeout=TIMEOUT) as ser, \
41.     open(OUT_FILENAME, 'w') as file:
42.     # output first line
43.     first_line = "Time/s\tCh1/nA\tCh2/nA\tCh3/nA\tCh4/nA\tCh5/nA\tCh6/nA\tT/°C"
44.     print(first_line)
45.     file.write(first_line + "\n")
46.     while(1): # Script will run until Ctrl+Z or Ctrl+C is pressed
47.         # Read data from serial port if available in the input buffer
48.         data = None
49.         if ser.inWaiting(): # is data in input buffer?
50.             data = ser.read()
```

Disclaimer: Not for medical, diagnostics and usage on humans. For evaluation use only.

```
51.         if data is not None:
52.             # operating with a ring array of size 25.
53.             # data is inserted at 0 and removed at 25. Therefore,
54.             # the data is inverted order in the array.
55.             # array ==> [newest_byte .... oldest_byte]
56.             # Each time the code parses the ring array is parsed until
57.             # a valid telegram is recognized.
58.             data_block.insert(0, data) # last byte in position 0
59.             nirvana = data_block.pop()
60.             del(nirvana)                # position 26 is deleted
61.             header = [b'\x04',          # Telegram header (inverted)
62.                       b'\x68',          # see documentation
63.                       b'\x13',
64.                       b'\x13',
65.                       b'\x68']
66.             cks = 0 # checksum
67.             # calculating checksum from byte 4 till second to last byte
68.             for x in [int.from_bytes(x, 'big') for x in data_block[2:-4]]:
69.                 cks = (cks + x) & 0xFF
70.             # validating header, end byte (x16) and checksum
71.             if data_block[-5:] == header \
72.                 and data_block[0] == b'\x16' \
73.                 and int.from_bytes(data_block[1], 'big') == cks:
74.                 # now inverting data train
75.                 # Useful data are from byte 5 until second to last byte.
76.                 data_inv = [x for x in data_block[2:-5]]
77.                 data_inv.reverse()
78.                 it = iter(data_inv) # iterator used to fetch 2 bytes
79.                 # next line turns 2 bytes into a 16-bit integer array
80.                 out_data = [
81.                     int.from_bytes(b''.join([x, next(it)]),
82.                                     byteorder='big',
83.                                     signed=True) for x in it]
84.                 # data to be stored
85.                 to_save = [str(x) for x in out_data] # 16-bit signed integers
86.                 # converting input data to currents in nanoamperes
87.                 # 50 nA ==> 32767 (2^15-1)
88.                 # applying gain and updating data as nanoamp.
89.                 gain = 50 / (2**15 - 1)
90.                 to_insert = [str(round(int(x) * gain, 3)) for x in to_save[0:6]]
91.                 # converting temperature to °C
92.                 temperature = str(round(float(to_save[6]) / 16, 3))
93.                 to_insert.append(temperature)
94.                 # generating timestamps
95.                 timestamp = round(time.time(), 4)
96.                 if start_timestamp is None:
97.                     start_timestamp = timestamp
98.                     delta_time = 0
99.                 else:
100.                    delta_time = timestamp - start_timestamp
101.                    to_insert.insert(0, str(round(delta_time, 1)))
102.                    # printing data on screen and output file
103.                    print("\t".join(to_insert))
104.                    file.write("\t".join(to_insert) + "\n")
105.                    # Next data point arrives in about 1.6 seconds
106.                    # So, resting for 1.4 seconds should be enough
107.                    time.sleep(1.4)
108.                    time.sleep(0.01) # Default polling time
```

Disclaimer: Not for medical, diagnostics and usage on humans. For evaluation use only.